

**\*\*Rapid Rhythm Reader does not detect heart attacks.**

**\*\*Rapid Rhythm Reader does not replace ordinary medical care.**

**\*\*The purpose of Rapid Rhythm Reader is to evaluate heart rhythm recordings obtained from mobile devices in persons without feelings of chest, arm, or jaw discomfort, shortness of breath, severe weakness, uncomfortable heart pounding, dizziness, or passing out. *If you have any of these feelings, you should seek medical attention at once.***

**\*\*Rapid Rhythm Reader is not a substitute for appropriate in-person or Telehealth care by a medical provider.**

## **How to send recordings for review:**

Although the *Apple Watch*® is used to obtain an ECG, the *iPhone*® is used to transmit an ECG recording for review.

1. On the paired iPhone, Select an ECG recording within the Apple Health app.
2. Tap the share icon.
3. Select Rapid Rhythm Reader.
4. Confirm in the dialog(s) that appear.

## **Where are my interpretations?**

On the home screen appears “My ECGs.” This is where a list of your ECG(s) that have been submitted will appear. An hourglass symbol indicates the recording has been submitted, but not yet reviewed. A red heart symbol indicates the ECG rhythm has been reviewed by the Cardiologist and is ready for your viewing.

## **I cannot see my interpreted ECGs:**

1. On the home screen entitled “My ECGs,” pull down the page until a refresh wheel appears.
2. A red heart symbol indicates an interpretation is ready for your viewing.
3. If an interpretation has not appeared within 24 hours of recording submission, email us at [help@healthalert.io](mailto:help@healthalert.io).

