

*****Rapid Rhythm Reader does not detect heart attacks.***

*****Rapid Rhythm Reader does not replace ordinary medical care.***

*****The purpose of Rapid Rhythm Reader is to evaluate heart rhythm recordings obtained from mobile devices in persons without feelings of chest, arm, or jaw discomfort, shortness of breath, severe weakness, uncomfortable heart pounding, dizziness, or passing out. If you have any of these feelings, you should seek medical attention at once.***

*****Rapid Rhythm Reader is not a substitute for appropriate in-person or Telehealth care by a medical provider.***

KardiaMobile for iOS-iPhone—Help Page

1—Follow instructions in KardiaMobile Quick Start Guide

2—Once ECG is completed, scroll down to the “more options” section and select download PDF

3—On password protection pop-up screen press “Skip.”

4—Press the share button in top right corner, and select the Rapid Rhythm Reader app.

5—Press “OK” on the sharing ECG pop-up screen.

6—Press “OK” on the next pop-up screen.